







**The Dolomite Centre Ltd.**  
 Unit 1, Anglian Business Park, Orchard Road,  
 Royston, Hertfordshire, SG8 5TW, UK

**T:** +44 (0)1763 242491  
**F:** +44 (0)1763 246125  
**E:** dolomitesales@syrris.com  
**W:** www.dolomite-centre.com

**The Dolomite Centre Limited**  
 Registered office: 27 Jarman Way,  
 Royston, Hertfordshire, SG8 5HW, UK  
 Company No. 05640084

## INFORMATION SHEET

Part name	Mitos Bottles	Part number	Various
-----------	---------------	-------------	---------

Name	Part Number	Description	Image
<b>Mitos Collection Bottle</b>	3000289	A 500ml glass bottle with an adaptor on the top for a Mitos Collection Needle (3000259). The lid is fitted with an o-ring to secure the Collection Needle and a vent-hole to prevent unwanted pressurisation.	
<b>Mitos Input Bottle 100ml</b>	3000290	A 100ml Duran bottle with connections for a Mitos Bottle Gas Hose and a 1/4"-28 thread. There is also a septum in the lid. The 1/4"-28 thread connects via a tube to a 5µm filter that is located at the base of the bottle. The maximum operating pressure is 1.25bar and the max temperature is 180°C.	
<b>Mitos Output Bottle 100ml</b>	3000293	A 100ml Duran bottle with a 1/4"-28 thread and a hole in the cap. Ideal for use as an output reservoir. The maximum operating pressure is 1.25bar and the maximum temperature is 180°C.	
<b>Mitos 250ml Multi-Port Bottle</b>	3000291	A 250ml Duran bottle with connections for a Mitos Bottle Gas Hose and four 1/4"-28 threads. There is also a septum in the lid. The four 1/4"-28 threads connect via a tube to a 5µm filter that is located at the base of the bottle. Maximum operating pressure is 1.25bar and the max temperature is 180°C.	
<b>Mitos Bottle Gas Hose</b>	3000292	This hose is ideal for connecting Mitos Input Bottle's to a gas supply.	